


Product Information	
Catalog No.	RA10938
Product name	InVivoMAb Anti-Human CD132/IL2RG (Iv0107)
Applications	ELISA, Neutralization
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	Iv0107
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	p64, Interleukin-2 receptor subunit gamma, gammaC, IL-2 receptor subunit gamma, CD132, IL-2R subunit gamma, IL-2RG, IL2RG, Cytokine receptor common subunit gamma
Accession	P31785
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human CD132/IL2RG (lv0107).

SDS-PAGE