


| Product Information       |                                                                                                                                                                                                                 |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.               | RA10963                                                                                                                                                                                                         |
| Product name              | Anti-DEX/Dextran Antibody (SAA0556)                                                                                                                                                                             |
| Applications              | ELISA                                                                                                                                                                                                           |
| Species reactivity        | General                                                                                                                                                                                                         |
| Host species              | Mouse                                                                                                                                                                                                           |
| Clonality                 | Monoclonal                                                                                                                                                                                                      |
| Isotype                   | IgG2a                                                                                                                                                                                                           |
| Clone ID                  | SAA0556                                                                                                                                                                                                         |
| Form                      | Liquid                                                                                                                                                                                                          |
| Storage buffer            | 0.01M PBS, pH 7.4.                                                                                                                                                                                              |
| Concentration             | 0.73 mg/ml                                                                                                                                                                                                      |
| Purity                    | >95% as determined by SDS-PAGE.                                                                                                                                                                                 |
| Purification              | Protein A/G purified from cell culture supernatant.                                                                                                                                                             |
| Target                    | (1, 6-alpha-D-Glucosyl)m, (1, 6-alpha-D-Glucosyl)m+1, (1, 6-alpha-D-Glucosyl)n, (1, 6-alpha-D-Glucosyl)n+1, 1, 6-alpha-D-Glucan, Dextran 40, Dextran 40, Dextran 70, Dextran 70, Dextran 75, dextrane, dextrans |
| Endotoxin level           | Please contact with the lab for this information.                                                                                                                                                               |
| Stability and Storage     | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.                                                    |
| Product Usage Information | ELISA: 1:2000-1:20000                                                                                                                                                                                           |
| Note                      | For research use only.                                                                                                                                                                                          |

### Image Data



SDS-PAGE for Anti-DEX/Dextran Antibody (SAA0556).

SDS-PAGE