


Product Information	
Catalog No.	RA10987
Product name	Anti-Vibrio cholerae LPS/Lipopolysaccharide Antibody (SAA0582)
Applications	ELISA
Species reactivity	Vibrio cholerae
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0582
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Lipopolysaccharide (LPS), O-specific polysaccharide (Vibrio cholerae serotype O1)
Accession	ChEBI: 16412
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Vibrio cholerae LPS/Lipopolysaccharide Antibody (SAA0582).

SDS-PAGE