


Product Information	
Catalog No.	RA11020
Product name	Anti-CH1&CH2/Allergen Fel d I-A&B Antibody (SAA0615)
Applications	ELISA, SPR
Species reactivity	Cat
Host species	Human
Clonality	Monoclonal
Isotype	IgG4
Clone ID	SAA0615
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CH2, AG4, Major allergen I polypeptide chain 2, Allergen Cat-1, Allergen Fel d I-B, Allergen FdI, Fel dI, Fel d 1-B, CH1, AG4, Major allergen I polypeptide chain 1, Allergen Cat-1, Allergen Fel d I-A, Allergen FdI, Fel dI, Fel d 1-A
Accession	P30438 & P30440
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-CH1&CH2/Allergen Fel d I-A&B Antibody (SAA0615).

SDS-PAGE