


Product Information	
Catalog No.	RA11059
Product name	Anti-mitF/aspF1 Antibody (SAA0654)
Applications	ELISA, WB
Species reactivity	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSCA1100) (Aspergillus fumigatus)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0654
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	mitF, aspF1, Ribonuclease mitogillin, Allergen Asp f I, Allergen I/a, IgE-binding ribotoxin, Major allergen Asp f 1, Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSCA1100) (Aspergillus fumigatus)
Accession	P67875
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-mitF/aspF1 Antibody (SAA0654).

SDS-PAGE