

Product Information	
Catalog No.	RA11061
Product name	Anti-Bet v 1/Allergen Bet v I-A Antibody (SAA0656)
Applications	ELISA
Species reactivity	Betula pendula (European white birch) (Betula verrucosa)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0656
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	BETVIA, BETVI, Major pollen allergen Bet v 1-A, Allergen Bet v 1-A, Betula pendula (European white birch) (Betula verrucosa)
Accession	P15494
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data

Anti-Bet v 1/Allergen Bet v I-A Antibody (SAA0656)
binds with Bet v 1/Allergen Bet v I-A

Detects Recombinant European white birch Bet v 1/Allergen Bet v I-A Protein, N-His (Cat No.: PR095012) in indirect ELISA.

