


Product Information	
Catalog No.	RA11081
Product name	Anti-PHLPI/Allergen Phl p I Antibody (SAA0676)
Applications	ELISA
Species reactivity	Phleum pratense (Common timothy)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0676
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	PHLPI, Pollen allergen Phl p 1, Allergen Phl p I, Phleum pratense (Common timothy)
Accession	P43213
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-PHLPI/Allergen Phl p I Antibody (SAA0676)

SDS-PAGE