


Product Information	
Catalog No.	RA11099
Product name	Anti-PHLPVa/Allergen Phl p Va Antibody (SAA0694)
Applications	ELISA, SPR
Species reactivity	Phleum pratense (Common timothy)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0694
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	PHLPVa, Pollen allergen Phl p 5a, Allergen Phl p Va, Phleum pratense (Common timothy)
Accession	Q40962
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-PHLPVa/Allergen Phl p Va Antibody (SAA0694)

SDS-PAGE