


Product Information	
Catalog No.	RA11155
Product name	Anti-PRUDU6 Antibody (SAA0750)
Applications	ELISA, WB
Species reactivity	Prunus dulcis (Almond) (Amygdalus dulcis)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0750
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	PRUDU6, Almond allergen Pru du 6, Allergen Pru du VI, 11S globulin, Almond Major Protein, Amandin, AMP, Pru-1, Prunin, Prunus dulcis (Almond) (Amygdalus dulcis)
Accession	E3SH28, Q43607
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-PRUDU6 Antibody (SAA0750)

SDS-PAGE