



Product Information	
Catalog No.	RA11168
Product name	Anti-Human GPRC5D Antibody (SAA0764)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0764
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	GPRC5D, G-protein coupled receptor family C group 5 member D
Accession	Q9NZD1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data




Flow Cytometry

Flow-cytometry using anti-Human GPRC5D antibody. GPRC5D Transfected 293T cells were stained with an irrelevant antibody (Green Histogram) or an anti-Human GPRC5D antibody monoclonal antibody (Catalog # HP458013, Blue Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a APC- conjugated goat anti-mouse antibody and cells analysed on a NovoCyte Flow Cytometer.n



SDS-PAGE

SDS-PAGE for Anti-Human GPRC5D Antibody (SAA0764).



Flow CytoMetry

Flow-cytometry using anti-human GPRC5D antibody. Untransfected cells (blue Histogram) and Transfected cells (Yellow Histogram) were stained with an anti-human GPRC5D monoclonal antibody (Catalog: HP458013) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a Goat Anti-Mouse IgG H&L Polyclonal Antibody, FITC (abInScience: MF690414) and cells analysed on a NovoCyte Flow Cytometer.