


Product Information	
Catalog No.	RA11169
Product name	Anti-H2-K1-peptide (SIINFEKL) Complex Antibody (SAA0757)
Applications	ELISA
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA0757
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.7 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	H-2 class I histocompatibility antigen, K-B alpha chain, H-2K(B), H2-K1, H2-K, Ovalbumin, Allergen Gal d II, Egg albumin, Plakalbumin, Gal d 2, SERPINB14
Accession	P01901 & P01012
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-H2-K1-peptide (SIINFEKL) Complex Antibody (SAA0757)

SDS-PAGE