


Product Information	
Catalog No.	RA11185
Product name	Anti-MPXV A29L Antibody (P2E12)
Applications	ELISA
Species reactivity	Monkeypox virus (strain Zaire-96-I-16) (MPX)
Host species	Human
Clonality	Monoclonal
Isotype	VHH-hFc
Clone ID	P2E12
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	A29L
Accession	NP_536566.1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-MPXV A29L Antibody (P2E12)

SDS-PAGE