


| Product Information   |                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | RA11210                                                                                                                                                      |
| Product name          | InVivoMAb Anti-DENV-2 Envelope protein E/EDIII (Iv0124)                                                                                                      |
| Applications          | ELISA, Neutralization                                                                                                                                        |
| Species reactivity    | Dengue virus type 2 (DENV-2)                                                                                                                                 |
| Host species          | Human                                                                                                                                                        |
| Clonality             | Monoclonal                                                                                                                                                   |
| Isotype               | IgG1                                                                                                                                                         |
| Clone ID              | Iv0124                                                                                                                                                       |
| Form                  | Liquid                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                           |
| Concentration         | 1 mg/ml                                                                                                                                                      |
| Purity                | >95% as determined by SDS-PAGE.                                                                                                                              |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                          |
| Target                | E glycoprotein, Envelope protein E, EDIII                                                                                                                    |
| Accession             | W5RZ25, P18356                                                                                                                                               |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                            |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| Note                  | For research use only. Not suitable for clinical or therapeutic use.                                                                                         |

### Image Data



SDS-PAGE for InVivoMAb Anti-DENV-2 Envelope protein E/EDIII (Iv0124).

SDS-PAGE