


Product Information	
Catalog No.	RA11250
Product name	InVivoMAb Anti-HRSV F/F protein Antibody (3M3)
Applications	ELISA, Neutralization
Species reactivity	Human respiratory syncytial virus (HRSV)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	3M3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	F, Fusion glycoprotein F0, Protein F
Accession	P12568, P11209, P03420, P13843, O36634
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-HRSV F/F protein Antibody (3M3)

SDS-PAGE