


Product Information	
Catalog No.	RA11256
Product name	Anti-Vibrio cholerae higB-2 Nanobody (SAA0851)
Applications	ELISA
Species reactivity	Vibrio cholerae serotype O1 (strain ATCC 39315 / El Tor Inaba N16961)
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA0851
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	Toxin HigB-2, higB-2, VC_A0468
Accession	Q9KMA6
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Vibrio cholerae higB-2 Nanobody (SAA0851)

SDS-PAGE