


Product Information	
Catalog No.	RA11274
Product name	Anti-CDV Protein N/Nucleoprotein Antibody (SAA0844)
Applications	ELISA, IF, WB
Species reactivity	Canine distemper virus (strain Onderstepoort) (CDV), Canine morbillivirus
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA0844
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Nucleoprotein, Protein N, Nucleocapsid protein, N
Accession	P04865
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-CDV Protein N/Nucleoprotein Antibody (SAA0844).

SDS-PAGE