


Product Information	
Catalog No.	RA11479
Product name	Anti-Human DPP6 Nanobody (SAA1245)
Applications	ELISA, FCM, SPR
Species reactivity	Human
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA1245
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	DPP VI, Dipeptidyl aminopeptidase-related protein, Dipeptidyl peptidase 6, DPP6, Dipeptidyl aminopeptidase-like protein 6, Dipeptidyl peptidase IV-like protein, Dipeptidyl peptidase VI, DPPX
Accession	P42658
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human DPP6 Nanobody (SAA1245)

SDS-PAGE