


Product Information	
Catalog No.	RA11482
Product name	Anti-Human HTT Nanobody (SAA1032)
Applications	ELISA
Species reactivity	Human
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA1032
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	Huntingtin, HD protein, Huntington disease protein, HTT, HD, IT15
Accession	P42858
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human HTT Nanobody (SAA1032)

SDS-PAGE