

Product Information	
Catalog No.	RA11487
Product name	Anti-Human CXCL12/SDF-1 Nanobody (SAA1273)
Applications	ELISA, RIA
Species reactivity	Human
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA1273
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	Intercrine reduced in hepatomas, hSDF-1, Stromal cell-derived factor 1, SDF-1, Pre-B cell growth-stimulating factor, SDF1, SDF1A, C-X-C motif chemokine 12, SDF1B, CXCL12, IRH, PBSF, hIRH
Accession	P48061
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.