

Product Information	
Catalog No.	RA11521
Product name	Anti-Human CDH17/Cadherin-17 Nanobody (SAA1325)
Applications	FCM, WB
Species reactivity	Human
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA1325
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.45 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	CDH17, Intestinal peptide-associated transporter HPT-1, LI-cadherin, Liver-intestine cadherin, Cadherin-17
Accession	Q12864
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CDH17/Cadherin-17 Nanobody (SAA1325)

SDS-PAGE