

Product Information	
Catalog No.	RA11675
Product name	Anti-SARS-CoV-2 S Protein Nanobody (SAA1122)
Applications	ELISA
Species reactivity	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA1122
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S1, Spike protein S2, Spike protein S2', S, 2
Accession	P0DTC2
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-SARS-CoV-2 S Protein Nanobody (SAA1122)

SDS-PAGE