

Product Information	
Catalog No.	RA11724
Product name	Anti-Clostridium botulinum botA/BOTOX Nanobody (SAA0932)
Applications	ELISA
Species reactivity	Clostridium botulinum (strain Hall / ATCC 3502 / NCTC 13319 / Type A)
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA0932
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	Botulinum neurotoxin type A, BoNT/A, Bontoxilysin-A, BOTOX, Botulinum neurotoxin A light chain, LC, 3.4.24.69, Botulinum neurotoxin A heavy chain, HC, botA, bna, CBO0806, CLC_0862
Accession	P0DPI1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Clostridium botulinum botA/BOTOX Nanobody (SAA0932)

SDS-PAGE