



Product Information	
Catalog No.	RA11777
Product name	Anti-Bacillus anthracis PA/pagA/PA-83 Antibody (SAA1376)
Applications	ELISA, SPR
Species reactivity	Bacillus anthracis
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA1376
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.5 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Protective antigen, PA, Anthrax toxins translocating protein, PA-83, PA83, Protective antigen PA-20, PA-20, PA20, Protective antigen PA-63, PA-63, PA63, pagA, pag, pXO1-110, BXA0164, GBAA_pXO1_0164
Accession	P13423
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data


**Anti-Bacillus anthracis PA/pagA/PA-83
Antibody (SAA1376) binds with PA**

Detects Bacillus anthracis PA/pagA/PA-83 in indirect ELISAs.



SDS-PAGE for Anti-Bacillus anthracis PA/pagA/PA-83 Antibody (SAA1376).

SDS-PAGE



Western blot analysis was performed using anti-Bacillus anthracis PA/pagA/PA-83 monoclonal antibody at 1ug/mL on various samples.
Lane 1: recombinant Bacillus anthracis PA/pagA/PA-83 (Catalog No: JN898022)
Lane 2: negative control

Western blot