


Product Information	
Catalog No.	RA11790
Product name	Anti-Bovine Serum Albumin/BSA Antibody (SAA1406)
Applications	ELISA
Species reactivity	Bovine
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	SAA1406
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.34 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Albumin, BSA, Bos d 6, ALB
Accession	P02769
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Bovine Serum Albumin/BSA Antibody (SAA1406).

SDS-PAGE