



Product Information	
Catalog No.	RA11794
Product name	Anti-Human MICA/B Antibody (20C6)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	20C6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	MICA, MIC-A, MHC class I polypeptide-related sequence A, PERB11.1, MIC-B, MHC class I polypeptide-related sequence B, MICB, PERB11.2
Accession	Q29983, Q29980
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human MICA/B Antibody (20C6).

SDS-PAGE



Flow CytoMetry

Flow-cytometry using anti-human MICA antibody. Untransfected cells (blue Histogram) and Transfected cells (Yellow Histogram) were stained with an anti-human MICA monoclonal antibody (Catalog: HS894013) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a Goat Anti-Mouse IgG H&L Polyclonal Antibody, FITC (abInScience: MF690414) and cells analysed on a NovoCyte Flow Cytometer.