


Product Information	
Catalog No.	RA11808
Product name	InVivoMAb Anti-Human VSIG4/CRlg Antibody (Iv0141)
Applications	ELISA, Neutralization
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH-Fc
Clone ID	Iv0141
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.98 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Protein Z39Ig, CRlg, V-set and immunoglobulin domain-containing protein 4, VSIG4, Z39IG
Accession	Q9Y279
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human VSIG4/CRlg Antibody (Iv0141).

SDS-PAGE