


Product Information	
Catalog No.	RA11810
Product name	Anti-HBV S/L-HBsAg/L glycoprotein Antibody (5a19)
Applications	ELISA, FuncS, WB
Species reactivity	Hepatitis B virus (HBV)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	5a19
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.67 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Large envelope protein S, L glycoprotein, L-HBsAg, LHB, Large surface protein, Major surface protein, preS1 domain, Pre-S1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-HBV S/L-HBsAg/L glycoprotein Antibody (5a19).

SDS-PAGE