


Product Information	
Catalog No.	RA11820
Product name	Anti-Influenza A virus Neuraminidase/NA protein Antibody (SAA1424)
Applications	ELISA
Species reactivity	Influenza A virus (A/environment/Shanghai/S1439/2013(H7N9))
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA1424
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Neuraminidase, 3.2.1.18, NA
Accession	S5MF06
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Influenza A virus Neuraminidase/NA protein Antibody (SAA1424)

SDS-PAGE