


Product Information	
Catalog No.	RA11867
Product name	Anti-4R-HNE/4-HNE Antibody (R310)
Applications	ELISA, WB
Species reactivity	General
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	R310
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HNE, 4-HNE, 4-hydroxynon-2-enal, 4-Hydroxy-2-nonenal, 4-Hydroxynonenal
Accession	ChEBI: 32585, 29343-52-0
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-4R-HNE/4-HNE Antibody (R310)

SDS-PAGE