


Product Information	
Catalog No.	RA11913
Product name	Anti-Plasmodium falciparum MTIP Antibody (SAA1445)
Applications	ELISA
Species reactivity	Plasmodium falciparum (isolate 3D7)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA1445
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Myosin A tail domain interacting protein, MTIP
Accession	Q8I4W8
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Plasmodium falciparum MTIP Antibody (SAA1445)

SDS-PAGE