


Product Information	
Catalog No.	RA11915
Product name	Anti-Opisthorchis viverrini Paramyosin Antibody (SAA1447)
Applications	ELISA, IHC, IP, WB
Species reactivity	Opisthorchis viverrini
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA1447
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Myosin heavy chain, Paramyosin
Accession	A0A1S8WLE5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-*Opisthorchis viverrini* Paramyosin Antibody (SAA1447)

SDS-PAGE