


Product Information	
Catalog No.	RA11916
Product name	Anti-Escherichia coli ptsH/hpr Antibody (Jel42)
Applications	ELISA
Species reactivity	Escherichia coli (strain K12)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Jel42
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	ptsH, hpr, Phosphocarrier protein HPr, Histidine-containing protein
Accession	P0AA04
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Escherichia coli ptsH/hpr Antibody (Jel42)

SDS-PAGE