


Product Information	
Catalog No.	RA11923
Product name	Anti-Human DLG4/PSD95 Antibody (K28/43)
Applications	IHC, IP, WB
Species reactivity	Human, Mouse, Rat
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	K28/43
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	PSD-95, PSD95, Disks large homolog 4, Synapse-associated protein 90, DLG4, Postsynaptic density protein 95, SAP-90, SAP90
Accession	P78352
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.


Image Data

Detects Human DLG4/PSD95 in indirect ELISAs.

**Anti-Human DLG4/PSD95 Antibody (K28/43)
binds with Human PSD95**



Bioactivity



SDS-PAGE for Anti-Human DLG4/PSD95 Antibody (K28/43).

SDS-PAGE