


Product Information	
Catalog No.	RA11956
Product name	Anti-Mouse MPO/Myeloperoxidase Antibody (6D6)
Applications	ELISA, IF, SPR, WB
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	6D6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Mpo, Myeloperoxidase, Myeloperoxidase light chain, Myeloperoxidase heavy chain
Accession	P11247
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse MPO/Myeloperoxidase Antibody (6D6)

SDS-PAGE