


Product Information	
Catalog No.	RA11970
Product name	Anti-Clostridium botulinum botA/BoNT Antibody (VH15)
Applications	ELISA, WB
Species reactivity	Clostridium botulinum
Host species	Human
Clonality	Monoclonal
Isotype	VHH-Fc
Clone ID	VH15
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	botA, atx, bonT, Botulinum neurotoxin type A, BoNT/A, Bontoxilysin-A, BOTOX, Botulinum neurotoxin type A1
Accession	P0DP10
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Clostridium botulinum botA/BoNT Antibody (VH15).

SDS-PAGE