



Product Information	
Catalog No.	RA11995
Product name	Anti-Vaccinia A27 Antibody (8E3#)
Applications	ELISA, FCM, WB
Species reactivity	Vaccinia
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	8E3#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.65 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Protein A27, VACWR150, A27L
Accession	P11258
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Vaccinia A27 Antibody (8E3#).

SDS-PAGE



Western blot analysis was performed using anti-A29L monoclonal antibody at 1ug/mL on various samples.

Lane 1: recombinant monkeypox virus A29L (Catalog No: VK415011)

Lane 2: negative control

Western blot