


Product Information	
Catalog No.	RA11996
Product name	Anti-Vaccinia A33R Antibody (A27D7)
Applications	ELISA, Neutralization, WB, IHC, IF
Species reactivity	Vaccinia
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	A27D7
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.42 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Protein A33, VACWR156, A33R
Accession	P68617
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Vaccinia A33R Antibody (A27D7).

SDS-PAGE