



Product Information	
Catalog No.	RA12030
Product name	InVivoMAb Anti-Human SNCA Antibody (Iv0173)
Applications	ELISA, Neutralization
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	Iv0173
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Alpha-synuclein, PARK1, Non-A4 component of amyloid precursor, NACP, Non-A beta component of AD amyloid, SNCA
Accession	P37840
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human SNCA Antibody (lv0173).

SDS-PAGE



Western blot analysis was performed using anti-SNCA monoclonal antibody at 1µg/ml on various samples.

Lane 1: recombinant human SNCA (Catalog No: HW755012)

Lane 2: negative control

Western blot