

Product Information	
Catalog No.	RA12048
Product name	InVivoMAb Anti-Mouse CD279/PDCD1/PD1 Antibody (RMP1-14)
Applications	Neutralization
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa (Silence)
Clone ID	RMP1-14
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Programmed cell death protein 1, Protein PD-1, hPD-1, PD1, PDCD1, CD279
Accession	Q02242
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.