


Product Information	
Catalog No.	RA12058
Product name	InVivoMAb Anti-Human MUC5AC Antibody (Iv0182)
Applications	ELISA, IHC, Neutralization
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Iv0182
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Mucin-5 subtype AC, tracheobronchial, Gastric mucin, Mucin-5AC, TBM, MUC-5AC, Tracheobronchial mucin, Major airway glycoprotein, MUC5AC, MUC5
Accession	P98088
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human MUC5AC Antibody (Iv0182).

SDS-PAGE