


Product Information	
Catalog No.	RA12065
Product name	InVivoMAb Anti-Human β -tryptase/TPSAB1 Antibody (31A)
Applications	Blocking, ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	31A
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TPSB1, TPS2, Tryptase I, Tryptase-1, TPSAB1, Tryptase alpha-1, Tryptase alpha/beta-1, TPS1
Accession	Q15661
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.


Image Data

Detects Human β -tryptase/TPSAB1 in indirect ELISAs.

Anti-Human β -tryptase/TPSAB1 Antibody (31A)
binds with TPSAB1




Bioactivity



SDS-PAGE for InVivoMAb Anti-Human β -tryptase/TPSAB1 Antibody (31A).

SDS-PAGE



Western blot analysis was performed using anti-TPSAB1/Tryptase-1 monoclonal antibody at 1ug/mL on various samples.

Lane 1: recombinant human TPSAB1/Tryptase-1 (Catalog No: HS931011)

Lane 2: negative control

Western blot