


Product Information	
Catalog No.	RA12076
Product name	Anti-Human FABP4/aP2 Antibody (SAA1478)
Applications	Blocking, ELISA
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA1478
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Adipocyte-type fatty acid-binding protein, Fatty acid-binding protein 4, ALBP, Adipocyte lipid-binding protein, FABP4, A-FABP, Fatty acid-binding protein, adipocyte, AFABP
Accession	P15090
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human FABP4/aP2 Antibody (SAA1478)

SDS-PAGE