


Product Information	
Catalog No.	RA12082
Product name	InVivoMAb Anti-Human TGM2 Antibody (52F3)
Applications	Blocking, ELISA, FCM, SPR
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	52F3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Protein-glutamine deamidase TGM2, TGase II, Heart G alpha(h), Erythrocyte transglutaminase, Protein-glutamine gamma-glutamyltransferase 2, TG(C), Transglutaminase C, Protein-glutamine serotonyltransferase TGM2, tTG, tTGase, Transglutaminase-2, Protein-glutamine noradrenalinyltransferase TGM2, TGM2, TG2, TGase C, Protein G alpha(h), Protein-glutamine histaminyltransferase TGM2, Tissue transglutaminase, hhG alpha(h), G(h), hTG2, Transglutaminase H, Protein-glutamine dopaminyltransferase TGM2, TGC, Isopeptidase TGM2, TGase-2, TGase H, Transglutaminase II
Accession	P21980
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human TGM2 Antibody (52F3).

SDS-PAGE