

Product Information	
Catalog No.	RA12087
Product name	InVivoMAb Anti-HSV1/HHV1 gD/Envelope glycoprotein D (Fd138-80)
Applications	Blocking, ELISA, FCM
Species reactivity	Human herpesvirus 1 (strain F) (HHV-1) (Human herpes simplex virus 1)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	Fd138-80
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Envelope glycoprotein D, gD, gD, US6
Accession	Q05059
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.