


Product Information	
Catalog No.	RA12095
Product name	Anti-Human GNAI1 & GNB1 Antibody (mAb16)
Applications	ELISA, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	mAb16
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	GNAI1, Guanine nucleotide-binding protein G(i) subunit alpha-1, Adenylate cyclase-inhibiting G alpha protein, Gi1 protein alpha subunit, G protein subunit alpha i1, Gαi1GNB1, Guanine nucleotide-binding protein G(l)/G(S)/G(T) subunit beta-1, Transducin beta chain 1, G protein subunit beta 1, Gβ1
Accession	P63096 & P62873
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human GNAI1 & GNB1 Antibody (mAb16)

SDS-PAGE