

Product Information	
Catalog No.	RA12107
Product name	Human Rheumatoid Factor IgM Antibody (2A2)
Applications	Rheumatoid Factor (RF) Reference Standard
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgM
Clone ID	2A2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein L purified from cell culture supernatant.
Target	Rheumatoid Factor, RF, rheumatoid factors (RFs)
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.