


Product Information	
Catalog No.	RA12138
Product name	Anti-Human ROBO1 Antibody (SAA1494)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1494
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	H-Robo-1, DUTT1, Roundabout homolog 1, Deleted in U twenty twenty, ROBO1
Accession	Q9Y6N7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.


Image Data



Western blot

Western blot analysis was performed using anti-ROBO1 monoclonal antibody at 1ug/mL on various samples.
Lane 1: recombinant human ROBO1 (Catalog No: HP485011)
Lane 2: negative control

**Anti-Human ROBO1 Antibody (SAA1494)
binds with ROBO1**



Bioactivity

Detects recombinant ROBO1 (Catalog No: HP485011) in indirect ELISA.