


Product Information	
Catalog No.	RA12244
Product name	Anti-Human PTN Antibody (SAA1600)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1600
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Heparin-binding neurite outgrowth-promoting factor 1, Heparin-binding brain mitogen, PTN, HBGF-8, Pleiotrophin, Heparin-binding growth factor 8, HBBM, HBNF-1, NEGF1, Osteoblast-specific factor 1, HBNF, HBNF1, Heparin-binding growth-associated molecule, HB-GAM, OSF-1, Heparin-binding neurite outgrowth-promoting factor
Accession	P21246
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human PTN Antibody (SAA1600)

SDS-PAGE