


Product Information	
Catalog No.	RA12261
Product name	Anti-Human EPHB2 Antibody (SAA1617)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1617
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	EPHT3, EPTH3, ERK, Developmentally-regulated Eph-related tyrosine kinase, Tyrosine-protein kinase receptor EPH-3, hEK5, Ephrin type-B receptor 2, Tyrosine-protein kinase TYRO5, EK5, HEK5, EPHB2, ELK-related tyrosine kinase, TYRO5, Renal carcinoma antigen NY-REN-47, EPH tyrosine kinase 3, DRT, EPH-like kinase 5
Accession	P29323
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human EPHB2 Antibody (SAA1617)

SDS-PAGE