


Product Information	
Catalog No.	RA12350
Product name	Anti-Human ACVR1C Antibody (SAA1706)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1706
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Activin receptor-like kinase 7, ALK7, ACVR1C, Activin receptor type IC, ALK-7, Activin receptor type-1C, ACTR-IC
Accession	Q8NER5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human ACVR1C Antibody (SAA1706)

SDS-PAGE